

1. SYSTEM INITIALIZATION AND SETUP

ALGORITHM: Initialize_System

```
BEGIN
  // Check GPU availability
  device ← get_gpu_device_name()
  IF device ≠ 'GPU:0' THEN
    RAISE SystemError('GPU not found')
  END IF

  // Set random seeds for reproducibility
  SET tensorflow_seed = 666
  SET numpy_seed = 666

  // Initialize hyperparameters
  BATCH_SIZE ← 64
  EPOCHS ← 50
  IMAGE_SIZE ← (224, 224)
  NUM_CLASSES ← 25
  LEARNING_RATE ← 2e-5
END
```

2. DATA AUGMENTATION PIPELINE

ALGORITHM: Data_Augmentation

```
FUNCTION flip_random_crop(image):
  image ← apply_random_horizontal_flip(image)
  image ← random_crop(image, size=(224, 224, 3))
  RETURN image

FUNCTION color_jitter(image, strength=[0.4, 0.4, 0.4, 0.1]):
  image ← adjust_brightness(image, delta=0.8*strength[0])
  image ← adjust_contrast(image, lower=1-0.8*strength[1], upper=1+0.8*strength[1])
  image ← adjust_saturation(image, lower=1-0.8*strength[2], upper=1+0.8*strength[2])
  image ← adjust_hue(image, delta=0.2*strength[3])
  image ← clip_values(image, min=0, max=255)
  RETURN image

FUNCTION color_drop(image):
```

```
image ← convert_to_grayscale(image)
image ← tile_to_rgb(image)
RETURN image
```

```
FUNCTION random_apply(func, image, probability):
  IF random() < probability THEN
    RETURN func(image)
  ELSE
    RETURN image
  END IF
```

```
FUNCTION custom_augment(image, label):
  image ← flip_random_crop(image)
  image ← random_apply(color_jitter, image, p=0.8)
  image ← random_apply(color_drop, image, p=0.2)
  RETURN image, label
```

3. DATASET PREPARATION

ALGORITHM: Load_Dataset

```
BEGIN
  // Define 25 classes for gastrointestinal images
  classes ← [
    cecum, ileum, retroflex_rectum, hemorrhoids, polyps,
    ulcerative_colitis_grade_0-1, ulcerative_colitis_grade_1,
    ulcerative_colitis_grade_1-2, ulcerative_colitis_grade_2,
    ulcerative_colitis_grade_2-3, ulcerative_colitis_grade_3,
    bbps_0-1, bbps_2-3, impacted_stool,
    dyed_lifted_polyps, dyed_resection_margins,
    pylorus, retroflex_stomach, z_line,
    barretts, barretts_short_segment,
    esophagitis_a, esophagitis_b-d,
    colon_adenocarcinoma, colon_benign_tissue
  ]

  // Create dataset dictionary mapping class_id to file_paths
  dataset ← {}
  FOR i = 0 TO 24 DO
    dataset[i] ← get_file_paths(classes[i])
  END FOR
END
```

```
FUNCTION combineDataAndLabel(width, height, dataset):  
  data ← empty_list()  
  labels ← empty_list()  
  
  FOR EACH class_id IN dataset DO  
    FOR EACH image_path IN dataset[class_id] DO  
      image ← load_image(image_path, size=(width, height))  
      data.append(image)  
      labels.append(class_id)  
    END FOR  
  END FOR  
END FOR  
  
data ← convert_to_array(data)  
labels ← convert_to_array(labels)  
RETURN data, labels
```

4. MOBILENETV3 MODEL ARCHITECTURE

ALGORITHM: Build_MobileNetV3_Model

```
FUNCTION create_model(num_classes, input_shape):  
    // Load pre-trained MobileNetV3Large model  
    base_model ← MobileNetV3Large(  
        weights='imagenet',  
        include_top=False,  
        input_shape=input_shape,  
        pooling='avg'  
    )  
  
    // Freeze base model layers initially  
    base_model.trainable ← False  
  
    // Build complete model  
    model ← Sequential([  
        Input(shape=input_shape),  
        base_model,  
        BatchNormalization(),  
        Dropout(rate=0.2),  
        Dense(units=num_classes, activation='softmax')  
    ])  
  
    RETURN model
```

5. MODEL TRAINING PROCESS

ALGORITHM: Train_Model

```
BEGIN  
    // Phase 1: Transfer Learning (Frozen Base)  
    model ← create_model(num_classes, input_shape)  
  
    // Configure optimizer  
    optimizer ← Adam(learning_rate=2e-5)  
  
    // Compile model  
    model.compile(  
        optimizer=optimizer,  
        loss='categorical_crossentropy',  
        metrics=['accuracy']
```

```

)

// Setup callbacks
callbacks ← [
    ReduceLROnPlateau(factor=0.2, patience=5),
    EarlyStopping(patience=10, restore_best_weights=True),
    ModelCheckpoint('best_model.h5', save_best_only=True)
]

// Train initial model
FOR epoch = 1 TO 20 DO
    train_batch ← get_next_batch(train_data, batch_size)
    augmented_batch ← custom_augment(train_batch)
    loss ← model.train_on_batch(augmented_batch)
    validate_model(model, validation_data)
END FOR

// Phase 2: Fine-tuning (Unfreeze Top Layers)
unfreeze_layers ← last_100_layers(base_model)
FOR layer IN unfreeze_layers DO
    layer.trainable ← True
END FOR

// Continue training with fine-tuning
optimizer ← Adam(learning_rate=1e-5) // Lower learning rate
model.compile(optimizer, loss, metrics)

FOR epoch = 21 TO 50 DO
    train_model_epoch(model, train_data, validation_data)
END FOR
END

```

6. MODEL EVALUATION AND METRICS

ALGORITHM: Evaluate_Model

```
FUNCTION evaluate_performance(model, test_data, test_labels):  
    // Make predictions  
    predictions ← model.predict(test_data)  
    predicted_classes ← argmax(predictions, axis=1)  
  
    // Calculate metrics  
    accuracy ← calculate_accuracy(predicted_classes, test_labels)  
    precision ← calculate_precision(predicted_classes, test_labels, average='weighted')  
    recall ← calculate_recall(predicted_classes, test_labels, average='weighted')  
    f1_score ← calculate_f1_score(predicted_classes, test_labels, average='weighted')  
  
    // Generate confusion matrix  
    confusion_matrix ← compute_confusion_matrix(test_labels, predicted_classes)  
  
    // Generate per-class classification report  
    classification_report ← generate_classification_report(  
        test_labels, predicted_classes, class_names  
    )  
  
    RETURN {  
        'accuracy': accuracy,  
        'precision': precision,  
        'recall': recall,  
        'f1_score': f1_score,  
        'confusion_matrix': confusion_matrix,  
        'classification_report': classification_report  
    }
```

7. MODEL EXPLAINABILITY (GRAD-CAM)

ALGORITHM: Generate_Explainability_Visualizations

```
FUNCTION generate_gradcam_heatmap(model, image, class_id):  
    // Prepare model for Grad-CAM  
    grad_model ← create_gradient_model(model, last_conv_layer)  
  
    // Forward pass and compute gradients  
    WITH gradient_tape AS tape:  
        conv_outputs, predictions ← grad_model(image)
```

```

    class_channel ← predictions[:, class_id]

    // Compute gradients
    grads ← tape.gradient(class_channel, conv_outputs)

    // Pool gradients and compute weighted activation
    pooled_grads ← global_average_pooling(grads)
    heatmap ← sum(conv_outputs * pooled_grads, axis=-1)
    heatmap ← relu(heatmap) // Apply ReLU
    heatmap ← normalize(heatmap, min=0, max=1)

    RETURN heatmap

FUNCTION visualize_all_cams(model, test_image, true_label):
    // Generate multiple CAM visualizations
    gradcam ← generate_gradcam(model, test_image, true_label)
    gradcam_plus ← generate_gradcam_plusplus(model, test_image, true_label)
    scorecam ← generate_scorecam(model, test_image, true_label)
    layercam ← generate_layercam(model, test_image, true_label)

    // Create visualization grid
    display_grid([
        original_image,
        overlay(original_image, gradcam),
        overlay(original_image, gradcam_plus),
        overlay(original_image, scorecam),
        overlay(original_image, layercam)
    ])

```

8. MAIN EXECUTION PIPELINE

ALGORITHM: Main_Pipeline

BEGIN

// Step 1: Initialize system

Initialize_System()

// Step 2: Load and prepare dataset

dataset ← Load_Dataset()

data, labels ← combineDataAndLabel(224, 224, dataset)

// Step 3: Split data

x_train, x_temp, y_train, y_temp ← train_test_split(data, labels, test_size=0.2)

x_val, x_test, y_val, y_test ← train_test_split(x_temp, y_temp, test_size=0.5)

// Step 4: One-hot encode labels

y_train ← to_categorical(y_train, num_classes=25)

y_val ← to_categorical(y_val, num_classes=25)

y_test_categorical ← to_categorical(y_test, num_classes=25)

// Step 5: Create data pipelines

train_dataset ← create_tf_dataset(x_train, y_train)

train_dataset ← train_dataset.map(custom_augment).batch(64).prefetch()

val_dataset ← create_tf_dataset(x_val, y_val).batch(64)

// Step 6: Build and train model

model ← create_model(25, (224, 224, 3))

Train_Model(model, train_dataset, val_dataset)

// Step 7: Evaluate model

metrics ← evaluate_performance(model, x_test, y_test)

PRINT "Accuracy:", metrics.accuracy

PRINT "Precision:", metrics.precision

PRINT "Recall:", metrics.recall

PRINT "F1 Score:", metrics.f1_score

// Step 8: Generate visualizations

plot_confusion_matrix(metrics.confusion_matrix, class_names)

visualize_predictions(model, x_test, y_test, sample_count=16)

// Step 9: Generate explainability visualizations

FOR sample IN test_samples DO

 visualize_all_cams(model, sample.image, sample.label)


```
END FOR
```

```
// Step 10: Save final model
```

```
save_model(model, 'mobilenetv3_gastrointestinal_classifier.h5')
```

```
END
```

APPENDIX: KEY DATA STRUCTURES

STRUCTURE: Class_Mapping

Class_ID → Class_Name:

- 0 → cecum
- 1 → ileum
- 2 → retroflex_rectum
- 3 → hemorrhoids
- 4 → polyps
- 5 → ulcerative_colitis_grade_0-1
- 6 → ulcerative_colitis_grade_1
- 7 → ulcerative_colitis_grade_1-2
- 8 → ulcerative_colitis_grade_2
- 9 → ulcerative_colitis_grade_2-3
- 10 → ulcerative_colitis_grade_3
- 11 → bbps_0-1
- 12 → bbps_2-3
- 13 → impacted_stool
- 14 → dyed_lifted_polyps
- 15 → dyed_resection_margins
- 16 → pylorus
- 17 → retroflex_stomach
- 18 → z_line
- 19 → barretts
- 20 → barretts_short_segment
- 21 → esophagitis_a
- 22 → esophagitis_b-d
- 23 → colon_adenocarcinoma
- 24 → colon_benign_tissue

STRUCTURE: Model_Architecture

Input_Layer: (224, 224, 3)

Base_Model: MobileNetV3Large (pre-trained on ImageNet)

Global_Average_Pooling: Reduce spatial dimensions

Batch_Normalization: Normalize activations

Dropout: Rate=0.2 for regularization

Dense_Output: 25 units with softmax activation

STRUCTURE: Training_Configuration

Batch_Size: 64

Initial_Epochs: 20 (frozen base)

Fine_Tuning_Epochs: 30 (unfrozen top layers)

Initial_Learning_Rate: 2e-5

Fine_Tuning_Learning_Rate: 1e-5

Optimizer: Adam

Loss_Function: Categorical_Crossentropy

Metrics: ['accuracy']

Callbacks: [ReduceLROnPlateau, EarlyStopping, ModelCheckpoint]